

Lab-Tips® - LTP125 Large Knitted Polyester Swab

Cleanroom laundered double layer 100% polyester knit swab

The LTP125 large knitted polyester swab is an for cleaning or sampling broad flat surfaces. Its excellent choice for general purpose cleaning or cleanroom laundered knitted head provide ultra surface sampling. Its durable flat head and rigid -low NVR and ion levels. internal paddle provide maximum surface area



- Cleanroom laundered 100% continuous filament polyester knit head
- Long polypropylene handle with a rigid paddle head
- Thermally bonded head
- Laundered and packaged in Class 10 (ISO Class 4) cleanroom

Benefits

- Critically low particles, fibers, ions and extractables
- High abrasion resistance while soft and non-abrasive
- Excellent chemical resistance
- No contaminating adhesives and good sorbency
- Autoclavable
- Free from Silicone, Amides and DOP

Environmental

Complies with RoHS and REACH requirements

Applications

- Designed for the highest level of contamination control in critical processing applications
- Cleaning of sensitive equipment such as print heads
- Sampling surfaces for cleaning validation
- Applying lubricants and other liquids
- Removing excess materials
- Cleaning with solvents

Physical Characteristics

Head Material	Knitted polyester Handle Thickness		2.8mm {0.110"}		
Head Width	12.7mm	{0.50"}	Handle Length	100mm {3.94"}	
Head Thickness	5mm	{0.196"}	Total Swab Length	125mm {4.92"}	
Head Length	25mm	{0.984"}	Head Bond	Thermal	
Handle Material	Polypropylene		Handle Text	Berkshire	
Handle Width	5.5mm	{0.216"}			

Other Knitted Polyester tip swabs

- LTP1465 Long handle
- LTP 70F Short handle flexible head
- LTP70R Short handle rigid head

www.berkshire.com

Contact: Tel 1 800 242 7000 info@berkshire.com

Americas	Tel 1 413 528 2602	Fax 1 413 528 2614	info@berkshire.com
Europe	Tel 44 (0) 870 757 2877	Fax 44 (0) 870 757 2878	enquiries@berkshire.uk.com
SE Asia	Tel 65 6252 4313	Fax 65 6252 4312	enquiries@berkshire.com.sg
Japan	Tel 81 3 5827 2380	Fax 81 35827 2382	master@berkshire.co.jp



Technical Data:

ATTRIBUTE		LTP125	TEST METHOD
Non-Volatile Residue	IPA mg/swab	0.075	IEST-RP.CC004.3, Sec 7.1.2
Ions (μg/swab)	Chloride	0.11	IEST-RP.CC004.3, Sec 7.2.2
	Sulfate	0.061	
	Nitrate	ND*	
	Phosphate	ND*	
	Fluoride	ND*	
	Potassium	0.010	
	Calcium	0.082	
	Sodium	0.032	
	Magnesium	0.011	

Notes:

- Data shown are typical values
- Due to differences in test methods applied and equipment utilized by different swab manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- * ND* = below the detection limit of the equipment.

Order Information:

Product	Item Number	Size	Swabs/pk	Pks/cs	Case size
Large Knitted Polyester	LTP125.5	125mm	100	5	220 x 160 X 105mm

Other Berkshire Products



Wipers



Mop Systems



Face Masks



Glove Liners



Documentation Systems



Swabs